

Glioma Cell Culture

Materials

MHM

name	dose	ratio
10X Hormone Mix	50 mL	1:10
L-glu	5 mL	1:100
p/s	5 mL	1:100
10X DMEM	50 mL	1:10
30% glucose	10 mL	1:50
7.5% NaHCO ₃	7.5 mL	1:66.7
HEPES	2.5 mL	1:200
ddH ₂ O	375 mL	-

hg medium

name	dose	ratio
MHM with p/s	50 mL	-
B27	1000 uL	1:50
FGF	50 uL	1:1000
EGF	50 uL	1:1000

Methods

Culture (from stock)

- defrost GBM stock in 37°C water bath
- MHM wash x1
 - centrifuge: 1200 rpm x5min
 - aspiration
- resuspend with 10mL hg medium, apply all to Ultra-low attachment T25
- after confirmation with microscope, transfer to incubator
- culture for 5 days (approx.)

Immunofluorescence

day1

- prepare **adherent 24-well dish for observation**
- aspirate fibro. in the wells
- apply hg medium 500uL each
- transfer culturing suspension 1mL each from T75 into wells
- 37°C incubation for 1hr

- confirm with microscope that nanospheres firmly adhere to the floor
- aspiration
- process with PFA x45min (approx.) in RT
- PBS wash x4
- process with 0.3% TX x5min (approx.) in RT
- PBS wash x2
- apply Primary Ab Mix
 - e.g.) Primary Ab Mix (for 1 well)

name	dose	ratio
Blocking One	500 uL	-
anti-xxx Mouse mAb	0.5 uL	1:1000
anti-yyy Rabbit mAb	0.5 uL	1:1000

- 4°C o/n with shading

day2

- PBS wash x3
- apply Secondary Ab Mix
 - e.g.) Secondary Ab Mix (for 1 well)

name	dose	ratio
Blocking One	500 uL	-
anti-Mouse Gt	0.5 uL	1:1000
anti-Rabbit Dk	0.5 uL	1:1000

- shade in RT x1hr (approx.)
- PBS wash x2
- apply 1:1000 hocheist

name	dose	ratio
PBS	500 uL	-
hocheist	0.5 uL	1:1000

- shade in RT x5min (approx.)
- PBS wash x2
- wait until wells get dry
- observe with KEYENCE